


# Imaging photodesorption from low temperature O<sub>2</sub> ice

Z. Sun,<sup>1,2</sup> R. A. Scheidsbach,<sup>1</sup> H. M. Cuppen,<sup>3</sup> H. V. J. Linnarts,<sup>4</sup> D. H. Parker<sup>1\*</sup>

1. Molecular and Laser Physics department, IMM, Radboud University, the Netherlands
2. School of Physics and Electronic Information, Anhui Normal University, China
3. Theoretical Chemistry department, IMM, Radboud University, the Netherlands
4. Laboratory for Astrophysics, Leiden Observatory, Leiden University, the Netherlands




## Abstract

A new mass spectrometry method, Atom Fragment Imaging (AFI-MS) provides for the first time both kinetic and internal energy data on mass selected photodesorption O atom products from molecular O<sub>2</sub> ice at 20K. Relative O : O<sub>2</sub> desorption yields found here are consistent with previous studies with the added surprise that a large fraction of the O<sub>2</sub> product is in the excited electronic states (a <sup>1</sup>Δ<sub>g</sub>, A' <sup>3</sup>Δ<sub>u</sub>) instead of the X<sup>3</sup>Σ<sub>g</sub><sup>-</sup> ground state. In addition, ozone is ejected from surface with VUV desorption.


## Motivation

In cold interstellar space (T=10-100 K) there is too much gas than expected; it may originate from photodesorption of ice-covered grains. What are the reaction mechanism when cold ice surfaces are exposed to UV/VUV light? In previous studies, REMPI-TOFMS spectroscopy was commonly used, yet only one-dimensional information could be obtained. Here we show that by combining REMPI detection and VMI (Velocity Map Imaging) technique in ice surface study, a new mass spectrometry method, Atom Fragment Imaging (AFI-MS) is developed, delivering rich three-dimensional velocity information to reveal photochemical pathways. O<sub>2</sub> is recently revealed to be the relevant in ISM and O<sub>2</sub> ice is a good starting system.

## Experiment setup and Methods




VMI records the recoil pattern for any TOF  $\Delta t$ , showing different origins of O(<sup>3</sup>P) atoms




## Results and Discussion

Beautiful images, even from ice!




## Conclusion

1. Combination of REMPI-VMI with ice is good: high sensitivity, selectivity and information content of REMPI-VMI is proven for ice surfaces by giving nice images;
2. VMI can greatly enhance TOF-MS because imaging pattern helps identify clearly the molecular precursors of O-atom signals, i.e, to reveal photochemical pathways on surfaces;
3. With UV light irradiation, Two-Photon Absorption (TPA) is necessarily taking place; with VUV light, One-Photon Absorption (OPA) takes place. Based on the characterizing measurements of O<sub>2</sub>-Rare gas mixed ice, the TPA cross section of O<sub>2</sub> ice is estimated to be unexpectedly large;
4. With UV desorption, the primary process is O<sub>2</sub> photodissociation together with recombination of O atoms to produce energetic O atoms and O<sub>2</sub> metastable molecules, respectively; similarity between discharge in O<sub>2</sub> and photodesorption of O<sub>2</sub> ice suggests recombination is important;
5. TPA occurs below the top surface, thus ground-state O<sub>2</sub> is largely made due to "Kick-off" mechanism with UV desorption. However OPA is effective at the surface with VUV light, giving higher probability of recombination between O and O<sub>2</sub>, leading to ozone ejection. And geminate recombination of two O atoms is more likely to generate highly vibrational O<sub>2</sub> molecular products.



Two-Photon Absorption



One-Photon Absorption

## Acknowledgement:

We gratefully acknowledge NWO and NSFC for financial supports.